


PELLITec®

BLISTER PREVENTION PADS

WITH SILIPOS® TECHNOLOGY

92892
2 per pack

92893
10 per poly bag



FEATURES

- ✓ **ADHESIVE LAYER** - to adhere securely
- ✓ **SCUBA LAYER** - tear resistant and water resistant
- ✓ **MYLAR LAYER** - for strength and durability
- ✓ **GEL LAYER** - creates cushioning effect and relaxes friction
- ✓ **MYLAR LAYER** - for strength and durability
- ✓ **BALTEX LAYER** - promotes breathability

The unique and highly effective sandwich construction of PelliTec® Blister Prevention Pads serves to move with your foot and helps to reduce the blisters caused by friction. By sticking to the inside of your footwear and over the blister, the pad works to cushion and protect the affected area.

The Benefits of PelliTec®:

- Significantly reducing the occurrence of blisters by walkers, sports players and the military
- Relieving pressure and pain
- Wicking away moisture and reducing the risk of infection

For more information please contact **SILIPOS®**

7049 Williams Road Niagara Falls, NY 14304

Tel: (800) 229-4404 • (716) 283-0700

Fax: (716) 283-0600 • Email: info@silipos.com

WWW.SILIPOS.COM

