

# SILIPOS®

---

## MOBILITY

MADE IN THE  
**USA**

# COMFORTZONE™

## GEL LINERS


ComfortZone Gel Liners the highest level in quality, protection, and comfort. ComfortZone Gel Liners offer maximum protection with medical grade mineral oil gel. ComfortZone gel liner is designed to provide maximum mobility and protection. The gel liner is adaptable with most prosthetic socket designs and is washable and reusable.

- 3 ply, 3mm of exposed gel
- Prevents abrasion and skin breakdown
- Optimum cushion and protection
- Ease of mobility
- Not made with natural rubber latex and hypoallergenic
- BPA & Phthalate Free

 With CoolMax®  
fabric




Silipos.com  
1-800-229-4404  
7049 Williams Road  
Niagara Falls, NY 14304



### SINGLE SOCKET GEL LINER

- 3 ply with a single layer of CoolMax<sup>®</sup> fabric
- Tapered liner shape
- 3mm thick medical grade mineral oil gel (against skin)



### DOUBLE SOCKET GEL LINER

- 5 ply with two layers of CoolMax<sup>®</sup> fabric
- Tapered liner shape
- 3mm thick medical grade mineral oil gel (between fabric)

## Double / Single Socket Gel Liner Sizing Guidelines

ITEM # SINGLE	ITEM # DOUBLE	SIZE	GEL LENGTH SINGLE	GEL LENGTH DOUBLE	FABRIC LENGTH SINGLE & DOUBLE	TOE / TOP CIRCUMFERENCE
12725	12625	Short / Narrow	~10" (~25cm)	~10" (~25cm)	~12" (~30cm)	~9"-12" / 15"-18" (~23-30cm / 38-46cm)
12755	12655	Short / Regular	~10" (~25cm)	~10" (~25cm)	~12" (~30cm)	~12"-14" / 18"-22" (~30-36cm / 45-56cm)
12765	12665	Medium / Regular	~10" (~25cm)	~12" (~30cm)	~16" (~41cm)	~9"-12" / 15"-18" (~30-36cm / 45-56cm)