

SILIPOS[®]

MOBILITY

MADE IN THE
USA



DURAGEL[™]

SUSPENSION SLEEVE & LINERS

The best combination of comfort, protection, and durability. The Silipos liner system has improved durability and is designed for higher activity levels without compromising our world renowned mineral oil gels that provide your patients with extreme comfort and total skin care.

- Liners – Available in 3mm, 6mm, and 9mm gel thickness
- Liners – Tapered and Uniformed styles
- Sleeves – 3mm gel thickness
- Sleeves - Easy to Donn and Doff
- Nylon composite fabric provides extreme abrasion resistance and superior gel bond
- Gel formulation will not tear, delaminate, or thin out with stress, maintaining skin care and overall patient comfort
- Hypoallergenic and not made with natural rubber latex
- BPA & Phthalate Free

DuraGel[™] has been designed for patients from low activity with extended life to moderate and active patients in the K1, K2, K3, and K4 levels.



Silipos.com
1-800-229-4404
7049 Williams Road
Niagara Falls, NY 14304



DURAGEL BLACK WOOD SUSPENSION SLEEVE

18306 – Medium | 18307 – Medium+ | 18308 – Large | 18309 – Large+ | 18310 – Medium+ Distal+

DURAGEL CAMO SUSPENSION SLEEVE

18186 – Medium | 18187 – Medium+ | 18188 – Large | 18189 – Large+ | 18190 – Medium+ Distal+

DURAGEL REINFORCED SUSPENSION SLEEVES

18181 – Medium | 18182 – Medium+ | 18185 – Medium+ Distal+ | 18183 – Large | 18184 – Large+

DURAGEL PIN LINER

3MM Uniform - 18211 – Medium | 18212 – Medium+ | 18213 – Large | 18214 – Large+ | 18215 – XLarge
3MM Tapered - 18216 – Medium | 18217 – Medium+
6MM Uniform - 18201 – Medium | 18202 – Medium+ | 18203 – Large | 18204 – Large+ | 18205 – XLarge
6MM Tapered - 18206 – Medium | 18207 – Medium+ | 18208 – Large
9MM - 18221 – Medium | 18222 – Medium+ | 18223 – Large | 18224 – Large+

DURAGEL CUSHION LINER

3MM Uniform - 18251 – Medium | 18252 – Medium+ | 18253 – Large | 18254 – Large+ | 18255 – XLarge
3MM Tapered - 18256 – Medium | 18257 – Medium+
6MM Uniform - 18241 – Medium | 18242 – Medium+ | 18243 – Large | 18244 – Large+ | 18245 – XLarge
6MM Tapered - 18246 – Medium | 18247 – Medium+ | 18248 – Large
9MM - 18261 – Medium | 18262 – Medium+ | 18263 – Large | 18264 – Large+

Liner Sizing Guidelines*

*All sizes are approximate.

| SIZE | PROXIMAL CIRCUMFERENCE | | DISTAL CIRCUMFERENCE | |
|---------|------------------------|-----------------|----------------------|-----------------|
| | RELAXED | FULLY STRETCHED | RELAXED | FULLY STRETCHED |
| Medium | 9" (23cm) | 15" (38cm) | 8" (20cm) | 13" (33cm) |
| Medium+ | 13" (33cm) | 20" (51cm) | 9" (23cm) | 13" (33cm) |
| Large | 15" (38cm) | 21" (53cm) | 9" (23cm) | 16" (41cm) |
| Large+ | 16" (41cm) | 26" (66cm) | 12" (30cm) | 17" (42cm) |
| X-Large | 16" (41cm) | 29" (74cm) | 12" (30cm) | 21" (53cm) |