

SILIPOS®

MOBILITY

MADE IN THE
USA

SILOSHEATH™


ANTIBACTERIAL

Antibacterial SiloSheath is specially formulated with a tested and proven antimicrobial ingredient that kills up to 99.9% of bacteria and fungus*. The antibacterial ingredient is contained in the gel for built in odor protection. This product is washable and reusable.

- Controls Odor
- Kills up to 99.9% of bacteria and fungus*
- Contains Silver Ion
- Antibacterial ingredient is contained in the gel for complete protection
- 2 ply, 2mm gel
- Minimizes skin breakdown
- Reduces friction, shearing, and abrasion on residual limb
- Softens, soothes, and improves the appearance of skin
- Hypoallergenic and not made with natural rubber latex
- BPA & Phthalate Free

 With Silipos®
antimicrobial gel

Silipos.com
1-800-229-4404
7049 Williams Road
Niagara Falls, NY 14304



ANTIBACTERIAL SILOSHEATH

- 14415 - 1
- 14425 - 2
- 14435 - 3
- 14445 - 4
- 14455 - 5
- 14465 - 6
- 14475 - 7
- 14485 - 8



- Kills up to 99.9% of bacteria and fungus*
- Contains Silver Ion
- Antibacterial ingredient is contained in the gel for complete protection

Silosheath™ Antibacterial Sizing Guidelines

SIZE	DESCRIPTION	APPROX. GEL LENGTH	APPROX. LENGTH	APPROX. DISTAL CIRCUMFERENCE	APPROX. PROXIMAL CIRCUMFERENCE
1	Child	8" (20cm)	10"-13" (25-33cm)	6½"-9" (17-23cm)	6½"-10" (17-25cm)
2	Short / Narrow	10" (25cm)	15"-17" (38-43cm)	7"-10½" (18-27cm)	7"-12" (18-30.5cm)
3	Short / Regular	10" (25cm)	15"-17" (38-43cm)	8"-12½" (20-32cm)	8"-14" (20-35.5cm)
4	Medium / Narrow	13" (33cm)	17"-19" (43-48cm)	8"-10½" (20-32cm)	8"-14" (20-35.5cm)
5	Medium / Regular	10" (25cm)	17"-19" (43-48cm)	9"-14" (23-35.5cm)	9"-16" (23-40.5cm)
6	Large / Narrow	13" (33cm)	18"-21" (46-53cm)	9"-14" (23-35.5cm)	9"-16" (23-40.5cm)
7	Large / Regular	13" (33cm)	18"-21" (46-53cm)	11"-16" (28-41cm)	11"-18" (28-46cm)
8	Large / Wide	13" (33cm)	18"-21" (46-53cm)	12"-19" (30.5-48cm)	12"-21" (30.5-53cm)