

# SILIPOS®

---

## MOBILITY




# SILOSHEATH™

## ORIGINAL

SiloSheath provides optimal range of motion and Comfort. SiloSheath is made with Silipos medical grade mineral oil gel and nylon outer material. SiloSheath is available in Original, Double Cushion and Extra Life styles. The Extra Life Silosheath is two nylon sheaths combined, with the mineral oil gel lining in-between the two layers of fabric. This style is recommended for individuals who prefer their skin to have no direct contact with the gel and for individuals who experience a higher degree of friction between their residual limbs and the prosthesis. Additionally, the extra layer of nylon helps to make the sheath more durable. All SiloSheath products are washable and reusable with 2 ply, 2mm gel and the SiloSheath Double Cushion is 2 ply, 3mm gel.

- Minimizes skin breakdown
- Reduces friction, shearing, and abrasion on residual limb
- Cushions and protects
- Softens, soothes, and improves the appearance of skin
- Hypoallergenic and not made with natural rubber latex
- BPA & Phthalate Free



Silipos.com  
1-800-229-4404  
7049 Williams Road  
Niagara Falls, NY 14304

### ORIGINAL EXTRA LIFE

13315 SIZE 1  
13325 SIZE 2  
13335 SIZE 3  
13345 SIZE 4  
13355 SIZE 5  
13365 SIZE 6  
13375 SIZE 7  
13385 SIZE 8


### DOUBLE CUSHION

13415 SIZE 1  
13425 SIZE 2  
13435 SIZE 3  
13445 SIZE 4  
13455 SIZE 5  
13465 SIZE 6  
13475 SIZE 7  
13485 SIZE 8

### EXTRA LIFE

13815 SIZE 1  
13825 SIZE 2  
13835 SIZE 3  
13845 SIZE 4  
13855 SIZE 5  
13865 SIZE 6  
13875 SIZE 7  
13885 SIZE 8

### SiloSheath<sup>™</sup> Original Sizing Guidelines



SIZE	DESCRIPTION	APPROX. GEL LENGTH	APPROX. LENGTH	APPROX. DISTAL CIRCUMFERENCE	APPROX. PROXIMAL CIRCUMFERENCE
1	Child	8" (20cm)	10"-13" (25-34cm)	6½"-9½" (17-24cm)	6½"-10½" (17-26cm)
2	Short / Narrow	10" (25cm)	15"-17" (38-44cm)	7"-11" (18-28cm)	7"-12½" (18-32cm)
3	Short / Regular	10" (25cm)	15"-17" (38-44cm)	8"-13" (20-34cm)	8"-15" (20-38cm)
4	Medium / Narrow	13" (33cm)	17"-19" (43-49cm)	8"-13" (20-34cm)	8"-15" (20-38cm)
5	Medium / Regular	10" (25cm)	17"-19" (43-49cm)	9"-15" (23-37cm)	9"-17" (23-43cm)
6	Large / Narrow	13" (33cm)	18"-21" (46-54cm)	9"-15" (23-37cm)	9"-17" (23-43cm)
7	Large / Regular	13" (33cm)	18"-21" (46-54cm)	11"-17" (28-42cm)	11"-19" (28-49cm)
8	Large / Wide	13" (33cm)	18"-21" (46-54cm)	12"-20" (31-50cm)	12"-22" (31-56cm)